

Product datasheet for **MR221565L4V**

Lgmn (NM_011175) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lgmn (NM_011175) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lgmn
Synonyms:	A; AEP; AI746452; AU022324; Pr; Prsc1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011175
ORF Size:	1308 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221565).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011175.2
RefSeq Size:	1874 bp
RefSeq ORF:	1308 bp
Locus ID:	19141
UniProt ID:	O89017
Cytogenetics:	12 E



[View online »](#)

Gene Summary:

This gene encodes a member of the cysteine peptidase family C13 that plays an important role in the endosome/lysosomal degradation system. The encoded inactive preproprotein undergoes autocatalytic removal of the C-terminal inhibitory propeptide to generate the active endopeptidase that cleaves protein substrates on the C-terminal side of asparagine residues. Mice lacking the encoded protein exhibit defects in the lysosomal processing of proteins resulting in their accumulation in the lysosomes, and develop symptoms resembling hemophagocytic lymphohistiocytosis. [provided by RefSeq, Aug 2016]