

Product datasheet for **MR221616L4V**

Shisa9 (NM_028277) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Shisa9 (NM_028277) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Shisa9
Synonyms:	2700045P11Rik; Ckamp44
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_028277
ORF Size:	1272 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221616).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_028277.2 , NP_082553.2
RefSeq Size:	4844 bp
RefSeq ORF:	1275 bp
Locus ID:	72555
UniProt ID:	Q9CZN4
Cytogenetics:	16
Gene Summary:	Regulator of short-term neuronal synaptic plasticity in the dentate gyrus. Associates with AMPA receptors (ionotropic glutamate receptors) in synaptic spines and promotes AMPA receptor desensitization at excitatory synapses.[UniProtKB/Swiss-Prot Function]



[View online »](#)