

Product datasheet for **MR221655L3V**

Dyrk1a (NM_007890) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dyrk1a (NM_007890) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dyrk1a
Synonyms:	2310043O08Rik; D16Ertd272e; D16Ertd493e; Dyrk; ENSMUSG00000074897; Gm10783; mmb; Mnbh; Mp86
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_007890
ORF Size:	2289 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221655).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007890.2 , NP_031916.1
RefSeq Size:	6035 bp
RefSeq ORF:	2292 bp
Locus ID:	13548
UniProt ID:	Q61214
Cytogenetics:	16 55.3 cM



[View online »](#)

Gene Summary:

Dual-specificity kinase which possesses both serine/threonine and tyrosine kinase activities. May play a role in a signaling pathway regulating nuclear functions of cell proliferation. Modulates alternative splicing by phosphorylating the splice factor SRSF6 (By similarity). Exhibits a substrate preference for proline at position P+1 and arginine at position P-3. Has pro-survival function and negatively regulates the apoptotic process. Promotes cell survival upon genotoxic stress through phosphorylation of SIRT1. This in turn inhibits TP53 activity and apoptosis (PubMed:20167603).[UniProtKB/Swiss-Prot Function]