

Product datasheet for **MR221660L3V**

Taf4b (NM_001100449) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Taf4b (NM_001100449) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Taf4b
Synonyms:	105kDa; 2610524B04Rik; 4932409F03Rik; AW987595; Taf2c2; TAFII-105; TAFII105
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001100449
ORF Size:	2565 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221660).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001100449.1
RefSeq Size:	5149 bp
RefSeq ORF:	2568 bp
Locus ID:	72504
UniProt ID:	G5E8Z2
Cytogenetics:	18 A1



[View online »](#)

Gene Summary:

Cell type-specific subunit of the general transcription factor TFIID that may function as a gene-selective coactivator in certain cells. TFIID is a multimeric protein complex that plays a central role in mediating promoter responses to various activators and repressors. TAF4B is a transcriptional coactivator of the p53/RELA NF-kappa-B subunit. Involved in the activation of a subset of antiapoptotic genes including TNFAIP3. Through interaction with OCBA/POU2AF1, acts as a coactivator of B-cell-specific transcription. Plays a role in spermiogenesis and oogenesis.[UniProtKB/Swiss-Prot Function]