

## Product datasheet for **MR221679L3V**

### Slc16a3 (NM\_001038653) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc16a3 (NM_001038653) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slc16a3
Synonyms:	Mct3; Mct4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001038653
ORF Size:	1410 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221679).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001038653.1</a> , <a href="#">NP_001033742.1</a>
RefSeq Size:	2481 bp
RefSeq ORF:	1413 bp
Locus ID:	80879
UniProt ID:	<a href="#">P57787</a>
Cytogenetics:	11 E2



[View online »](#)

**Gene Summary:**

Proton-linked monocarboxylate transporter. Catalyzes the rapid transport across the plasma membrane of many monocarboxylates such as lactate, pyruvate, branched-chain oxo acids derived from leucine, valine and isoleucine, and the ketone bodies acetoacetate, beta-hydroxybutyrate and acetate (By similarity).[UniProtKB/Swiss-Prot Function]