

## Product datasheet for **MR221687L3V**

### **Bicd2 (NM\_029791) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Bicd2 (NM_029791) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Bicd2
Synonyms:	0610027D24Rik; 1110005D12Rik; AA408834; mKIAA0699
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_029791
ORF Size:	2460 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221687).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_029791.3</a>
RefSeq Size:	4634 bp
RefSeq ORF:	2463 bp
Locus ID:	76895
UniProt ID:	<a href="#">Q921C5</a>
Cytogenetics:	13 A5



[View online »](#)

**Gene Summary:**

Acts as an adapter protein linking the dynein motor complex to various cargos and converts dynein from a non-processive to a highly processive motor in the presence of dynactin. Facilitates and stabilizes the interaction between dynein and dynactin and activates dynein processivity (the ability to move along a microtubule for a long distance without falling off the track) (PubMed:11483508, PubMed:25035494, PubMed:24986880, PubMed:22956769). Facilitates the binding of RAB6A to the Golgi by stabilizing its GTP-bound form (PubMed:25962623). Regulates coat complex coatamer protein I (COPI)-independent Golgi-endoplasmic reticulum transport via its interaction with RAB6A and recruitment of the dynein-dynactin motor complex (PubMed:12447383, PubMed:25962623). Contributes to nuclear and centrosomal positioning prior to mitotic entry through regulation of both dynein and kinesin-1. During G2 phase of the cell cycle, associates with RANBP2 at the nuclear pores and recruits dynein and dynactin to the nuclear envelope to ensure proper positioning of the nucleus relative to centrosomes prior to the onset of mitosis (PubMed:20386726). [UniProtKB/Swiss-Prot Function]