

## Product datasheet for **MR221700L3V**

### Zpbp (NM\_001185153) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Zpbp (NM_001185153) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Zpbp
Synonyms:	4930486K01Rik; lam38; Sp38
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001185153
ORF Size:	1053 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221700).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001185153.1</a> , <a href="#">NP_001172082.1</a>
RefSeq Size:	3674 bp
RefSeq ORF:	1056 bp
Locus ID:	53604
Cytogenetics:	11 A1
Gene Summary:	Plays a role in acrosome compaction and sperm morphogenesis. Is implicated in sperm-oocyte interaction during fertilization.[UniProtKB/Swiss-Prot Function]



[View online »](#)