

Product datasheet for **MR221717L4V**

Nup35 (NM_027091) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Nup35
Synonyms:	35kDa; 2310006I24Rik; 5330402E05Rik; MP44; NO44
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_027091
ORF Size:	975 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR221717).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_027091.4 , NP_081367.1
RefSeq Size:	2732 bp
RefSeq ORF:	978 bp
Locus ID:	69482
UniProt ID:	Q8R4R6
Cytogenetics:	2 C3



Gene Summary:

Functions as a component of the nuclear pore complex (NPC). NPC components, collectively referred to as nucleoporins (NUPs), can play the role of both NPC structural components and of docking or interaction partners for transiently associated nuclear transport factors. May play a role in the association of MAD1 with the NPC (By similarity).[UniProtKB/Swiss-Prot Function]