

Product datasheet for MR221723L4

Clec4b1 (NM_001190310) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Clec4b1 (NM_001190310) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Clec4b1
Synonyms:	1810046I24Rik; Aplra2; Clec4b; DCAR; Dcar2; DCARbeta
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221723).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

Locus ID: 69810

Cytogenetics: 6 F2