

Product datasheet for **MR221783**

Pkib (NM_001190402) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Pkib (NM_001190402) Mouse Tagged ORF Clone

Tag: Myc-DDK

Symbol: Pkib

Synonyms: Pkibeta; Prkacn2

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >MR221783 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGGACAGATTCATCAGAGATGGCTGATGTGGAATCTGTGATCACCAGCTTTGCGTCCTCAGCAAGGG
CAGGCCGCCGAATGCCTTACCCGACATCCAGAGTTCACTGGCTACAAGTGGATCCTCTGATCTTCCACT
GAAACTGGAAGCGTTGGCCGTGAAGGAAGATGCAAAAACGAAGAATGAAGAGAAAGACCAAGGCCAACCA
AAAACACCCTAAATGAAGGAAAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR221783 protein sequence
Red=Cloning site Green=Tags(s)

MRTDSSEMADVESVITSFASSARAGRRLPDIQSSLATSGSSDLPLKLEALAVKEDAKTKNEEKDQGQP
KTPLNEGK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_001190402

ORF Size: 237 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001190402.1](#), [NP_001177331.1](#)

RefSeq Size: 5485 bp


RefSeq ORF: 237 bp

Locus ID: 18768

Cytogenetics: 10 29.22 cM

MW: 8.3 kDa

Product images:



Circular map for MR221783