

Product datasheet for **MR221789L4V**

Gas2l1 (NM_144560) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gas2l1 (NM_144560) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gas2l1
Synonyms:	4930500E24Rik; D0Jmb1; Gar22; TU-71.1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_144560
ORF Size:	2037 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221789).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_144560.3 , NP_653146.1
RefSeq Size:	2858 bp
RefSeq ORF:	2037 bp
Locus ID:	78926
UniProt ID:	Q8JZP9
Cytogenetics:	11 A1



[View online »](#)

Gene Summary:

This gene encodes a member of the Gas2 family, actin-associated proteins expressed at high levels in growth-arrested cells. The gene expression is negatively regulated by serum and growth factors. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2010]