

Product datasheet for **MR221803L3V**

Caly (NM_001190385) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Caly (NM_001190385) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Caly
Synonyms:	0710001P07Rik; 1110004A22Rik; Calcyon; Drd1ip
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001190385
ORF Size:	678 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221803).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001190385.1 , NP_001177314.1
RefSeq Size:	1260 bp
RefSeq ORF:	681 bp
Locus ID:	68566
UniProt ID:	Q9DCA7
Cytogenetics:	7 F4
Gene Summary:	Interacts with clathrin light chain A and stimulates clathrin self-assembly and clathrin-mediated endocytosis.[UniProtKB/Swiss-Prot Function]



[View online »](#)