

## Product datasheet for **MR221833L3V**

### Arhgef1 (NM\_008488) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Arhgef1 (NM_008488) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Arhgef1
Synonyms:	Lbcl2; Lsc
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008488
ORF Size:	2757 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221833).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_008488.2</a> , <a href="#">NP_032514.1</a>
RefSeq Size:	3270 bp
RefSeq ORF:	2760 bp
Locus ID:	16801
UniProt ID:	<a href="#">Q61210</a>
Cytogenetics:	7 A3



[View online »](#)

**Gene Summary:**

Seems to play a role in the regulation of RhoA GTPase by guanine nucleotide-binding alpha-12 (GNA12) and alpha-13 (GNA13) subunits. Acts as GTPase-activating protein (GAP) for GNA12 and GNA13, and as guanine nucleotide exchange factor (GEF) for RhoA GTPase. Activated G alpha 13/GNA13 stimulates the RhoGEF activity through interaction with the RGS-like domain. This GEF activity is inhibited by binding to activated GNA12. Mediates angiotensin-2-induced RhoA activation. Isoform 3 and isoform 4 do not homooligomerize and show an enhanced RhoGEF activity.[UniProtKB/Swiss-Prot Function]