

Product datasheet for **MR221835L3V**

Mthfd1l (NM_172308) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mthfd1l (NM_172308) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mthfd1l
Synonyms:	2410004L15Rik; AI647056; Fthfsdc1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_172308
ORF Size:	2931 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221835).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172308.4 , NP_758512.3
RefSeq Size:	3731 bp
RefSeq ORF:	2934 bp
Locus ID:	270685
UniProt ID:	Q3V3R1
Cytogenetics:	10 A1



[View online »](#)

Gene Summary:

May provide the missing metabolic reaction required to link the mitochondria and the cytoplasm in the mammalian model of one-carbon folate metabolism in embryonic an transformed cells complementing thus the enzymatic activities of MTHFD2.[UniProtKB/Swiss-Prot Function]