

Product datasheet for **MR221848**

Cmah (NM_001111110) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cmah (NM_001111110) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Cmah
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR221848 representing NM_001111110
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGATGGACAGGAAACAGACAGCTGAGACCCTGCTGACCCTGCTCCTGCTGAAGTTGCCAACCTCAAGG
 AAGGGATCAATTTTTTCGAAATAAGACTACTGGGAAAGAGTACATTTTATACAAGGAGAAGGACCATCT
 AAAGGCATGCAAGAACCTCTGCAAGCACCAGGGAGGCCTGTTTCATGAAAGACATCGAGGATTTAGATGGA
 AGGTCCGTTAAATGCACAAAGCACAACCTGGAAGTTAGACGTGAGCACCATGAAATATATCAACCCTCCAG
 GGAGCTTCTGTCAAGACGAGCTCGTTATTGAAATGGATGAAAACAATGGGCTTTCCCTGGTAGAACTGAA
 CCCTCCTAACCCCTGGGACTCTGATCCCAGGTCTCCTGAAGAATTAGCTTTTGGGAAAGTACAGATAACA
 TATCTCACTCATGCCTGCATGGACCTCAAGTTGGGAGACAAGCGAATGGTATTTGACCCTTGTTAATTG
 GCCCTGCTTTGCCCGAGGATGGTGGTGTACATGAGCCTCCATCTGACTGGTTGGAGAGGCTGTGCAA
 AGCAGACCTCATTATATCAGCCACATGCACTCAGACCACCTGAGCTACCCTACCCTGAAGCAGCTTTCC
 CAGAGACGACCAGACATTCCTTTATGTTGGCGACACAGAAAGGCCTGTGTTTTGGAACCTGGATCAGA
 GTGGCGTCGGGTTAACTAACATCAACGTGGTTCATTTGGAATATGGCAACAGGTAGACAAAAGTCTGCG
 GTTCATGATCTTGATGGACGGCTTCATCCTGAGATGGACACATGCATTATCGTGGAGTACAAAAGTTCAT
 AAAATACTCAACACAGTGGACTGCACCAGACCCAAATGGGGGAAGGCTTCTGAGAAAAGTTGCTCTAATGA
 TGAGTGATTTTCGAGGAGGTGCATCAGGCTTCCAATGACTTTCAGTGGTGGAAAATTTACTGAGGAATG
 GAAAGCCCAGTTCATTAAGGCTGAAAGAAGAAAGCTTCTGAATTAACAAGCTCAGCTGGTGAAGGACCTG
 CAGCCCCGAATCTACTGTCGGTTGCTGGTACTTTGTGGAGTCTACCCATCTGACAAGTACATTAAGG
 AAACAAACACCAAAAATGACCCAAATCAGCTCAACAATCTTATCAGGAAAACTCTGACGTGGTGACATG
 GACCCACGACCTGGCGCTGTCCTCGACCTTGGCAGGATGCTGAAGGACCAACAGACAGCAAGGGCATT
 GTGGAGCCTCCAGAGGGGACAAAGATTTACAAGGATTCCTGGGACTTTGGCCCGTACCTGGAGATCTTGA
 ATTCTGCTGTGAGAGATGAAATCTTCTGTCATTATCCTGGATTAAGAGTACTTCACGTGGGCTGGATT
 TAAGAATTACAACCTGGTGGTGCAGGATGATTGAAACAGATGAAGATTTTCAGCCCTTTTCTGGAGGGTAC
 GACTATCTGGTGGACTTTCTAGATTTATCCTTTCCGAAAGAAAGACCCAGCCGGGAGCATCCTTATGAAG
 AAATCCATAGCCGGTGGATGTCATCAGGTACGTGGTGAAGAACGGCCTGCTGTGGGATGATCTGTATAT
 TGGATTCCAGACCCGATTGCTGCGGGACCCTGATATATACCATCATCTGTTTTGGAAATCATTTTCAGATA
 AAACCTCCTTAACACCACCAACTGGAAGTCGTTCTCAATGCACTGTGAT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR221848 representing NM_001111110
 Red=Cloning site Green=Tags(s)


MMDRKQTAETLLTLSPAIEVANLKEGINFFRNKTTGKEYILYKEKDHLKACKNLCKHQGGLFMKDIEDLDG
 RSVKCTKHNWKLVDSTMKYINPPGSFCQDELVIEMDENGLSLVELNPPNPWDS DPRSPEELAFGEVQIT
 YLTHACMDLKLGDKRMVFD PWLIGPAFARGWLLHEPPSDWLERLCKADLIYISHMHSDDL SYPTLKQLS
 QRRPDIPYVGDTERPVFVNLDQSGVGLTNINVVPFGIQQVDKSLRFMILMDGVHPMDTCIIIVEYKGH
 KILNTVDCTRPNGRLEPKVALMMSDFAGGASGFPMTFSGGKFTTEWKAQFIKAERRKLLNYKAQLVKDL
 QPRIYCPFAGYFVESHPDKYIKETNTKNDPNQLNLLIRKNSDVVTWTPRPGAVLDLGRMLKDPTDSKGI
 VEPPEGTKIYKDSWDFGPYLEILNSAVRDEIFCHSSWIKEYFTWAGFKNYNLVVRMIETDEDFSPFPGGY
 DYLVDFLDL SFPKERPSREHPYEEIHSRVDIVRYVVKNGLLWDDLYIGFQTRLLRDPDIYHHLFWNHFQI
 KLPLTPPNWKSFLMHCD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM_001111110

ORF Size:	1731 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001111110.2 , NP_001104580.1
RefSeq Size:	9727 bp
RefSeq ORF:	1734 bp
Locus ID:	12763
UniProt ID:	Q61419
Cytogenetics:	13 A3.1
MW:	67.4 kDa
Gene Summary:	Sialic acids are components of carbohydrate chains of glycoconjugates and are involved in cell-cell recognition and cell-pathogen interactions. Catalyzes the conversion of CMP-N-acetylneuraminic acid (CMP-Neu5Ac) into its hydroxylated derivative CMP-N-glycolylneuraminic acid (CMP-Neu5Gc), a sialic acid abundantly expressed at the surface of many cells.[UniProtKB/Swiss-Prot Function]