

Product datasheet for **MR221863L4V**

Rgs9bp (NM_145840) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rgs9bp (NM_145840) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rgs9bp
Synonyms:	R9AP; Rgs9-1bp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145840
ORF Size:	714 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221863).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145840.3 , NP_665839.1
RefSeq Size:	6604 bp
RefSeq ORF:	714 bp
Locus ID:	243923
UniProt ID:	Q148R9
Cytogenetics:	7 B2



[View online »](#)

Gene Summary:

Regulator of G protein-coupled receptor (GPCR) signaling in phototransduction. Participates in the recovery phase of visual transduction via its interaction with RGS9-1 isoform. Acts as a membrane-anchor that mediates the targeting of RGS9-1 to the photoreceptor outer segment, where phototransduction takes place. Enhances the ability of RGS9-1 to stimulate G protein GTPase activity, allowing the visual signal to be terminated on the physiologically time scale. It also controls the proteolytic stability of RGS9-1, probably by protecting it from degradation.[UniProtKB/Swiss-Prot Function]