

## Product datasheet for **MR221865L4V**

### Mrvi1 (Irag1) (NM\_010826) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Mrvi1 (Irag1) (NM_010826) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Irag1
Synonyms:	BB115629; Irag; Mrvi1; R; Ris1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_010826
ORF Size:	2763 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221865).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_010826.4</a> , <a href="#">NP_034956.2</a>
RefSeq Size:	5958 bp
RefSeq ORF:	2766 bp
Locus ID:	17540
Cytogenetics:	7 57.9 cM



[View online »](#)

**Gene Summary:**

This gene is a putative tumor suppressor gene that is frequently disrupted by mouse AIDS-related virus (MRV). The encoded protein participates in signaling by nitric oxide (NO) to inhibit intracellular calcium release and platelet aggregation in cardiovascular tissue. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Oct 2013]