

Product datasheet for **MR221872L4V**

MIph (NM_053015) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MIph (NM_053015) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	MIph
Synonyms:	2210418F23Rik; 5031433I09Rik; AW228792; D1Wsu84e; I(1)-3Rk; I1Rk3; In; Slac-2a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053015
ORF Size:	1770 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221872).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053015.3 , NP_443748.2
RefSeq Size:	4484 bp
RefSeq ORF:	1773 bp
Locus ID:	171531
UniProt ID:	Q91V27
Cytogenetics:	1 45.73 cM
Gene Summary:	Rab effector protein involved in melanosome transport. Serves as link between melanosome-bound RAB27A and the motor protein MYO5A.[UniProtKB/Swiss-Prot Function]



[View online »](#)