

Product datasheet for **MR221876L3V**

Ipcef1 (NM_001170800) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ipcef1 (NM_001170800) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ipcef1
Synonyms:	A130090K04Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001170800
ORF Size:	1218 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221876).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001170800.1 , NP_001164271.1
RefSeq Size:	5561 bp
RefSeq ORF:	1221 bp
Locus ID:	320495
UniProt ID:	Q5DU31
Cytogenetics:	10 A1
Gene Summary:	Enhances the promotion of guanine-nucleotide exchange by PSCD2 on ARF6 in a concentration-dependent manner.[UniProtKB/Swiss-Prot Function]



[View online »](#)