

Product datasheet for **MR221882**

Bcl11a (NM_016707) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Bcl11a (NM_016707) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Bcl11a
Synonyms:	2810047E18Rik; BCL-11A; Ctip1; D930021L15Rik; Evi9; Evi9a; Evi9b; Evi9c; mKIAA1809
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>MR221882 representing NM_016707
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTCTCGCCGCAAGCAAGGCAAAACCCAGCACTTAAGCAAACGGGAATTCGCCCCGAACCTCTTGAAG
 CCATTCTTACAGATGATGAACCAGACCATGGCCGTTGGGAGCTCCAGAAGGGGACCACGACCTTCTCAC
 CTGTGGGCACTGCCAGATGAATTTCCCACTGGGGGACATTCTTATTTTTATCGAGCAGAAACGAAACAA
 TGCAATGGCAGCCTCTGCTTAGAAAAAGGTGTGGATAAGCCGCCTTCCCCTTCTCCATCGAGATGAAAA
 AGGCATCCAATCCTGTGGAGTTGGCATCCAGGTACGCCAGAGGATGACGATTGTTTATCAACGTCATC
 TAGAGGAATTTGCCCCAAACAGGAACACATAGCAGATAAACTTCTGCACTGGAGGGGCTGTCTCTCT
 CGGTCTGCACACGGAGCTCTAATCCCCACGCCCGGGATGAGTGCAGAATATGCCCCGAGGGTATTTGTA
 AAGATGAGCCAGCAGCTACACATGTACAATTGCAAACAGCCATTACCAGTGCATGGTTTCTCTTGCA
 ACACGCACAGAACTCATGGATTAAGAATCTACTTAGAAAGTGAACACGGAAGTCCCCTGACCCCGCGG
 GTTGGTATCCCTTCAGGACTAGGTGCAGAATGTCTTCCAGCCACCTCTCCATGGGATTCATATTGCAG
 ACAATAACCCCTTTAACCTGCTAAGAATACCAGGATCAGTATCGAGAGAGGCTTCCGGCCTGGCAGAAGG
 GCGCTTTCCACCACTCCCCCTGTTAGTCCACCACCGAGACATCACTTGGACCCCCACCGCATAGAG
 CGCCTGGGGGCGGAAGAGATGGCCCTGGCCACCATCACCCGAGTGCCTTTGACAGGGTGTGCGTTGA
 ATCCAATGGCTATGGAGCCTCCCGCCATGGATTCTCTAGGAGACTTAGAGAGCTGGCAGGGAACACGTC
 TAGTCCACCGCTGTCCCAGGCCGCGCCAGTCTATGCAAAGTTACTGCAACCATTCAGCCAGGTAGC
 AAGCCACCCTTCTGGGACGCCCCCTCCCTCTCTGCAATCCGCCCTCTCCCTCCCAACCCCGG
 TCAAGTCCAAGTCATGCGAGTTCTGCGCAAGACGTTCAAATTTAGAGCAACTTGGTGGTTCACCGACG
 CAGCCATACTGGTGAGAAGCCCTATAAGTGCAACCTGTGCGACCACGCGTGCACACAGGCCAGCAAGCTG
 AAGCGTCACATGAAGACACATGCACAAATCGTCCCCCATGACAGTCAAGTCCGACGATGGCCTCTCCA
 CAGCCAGCTCCCCGAACCTGGTACCAGCGACCTGGTGGGACGCGCCAGCAGTGCCTCAAGTCAGTGGT
 GGCCAAGTTCAAGAGTGAGAACGACCCCACTTGATCCAGAGAACGGGGATGAGGAGGAAGAGGAGGAC
 GACGAGGAAGAAGAAGAGGAGGAAGAGGAGGAGGAGGAGCTGACGGAGAGCGAGAGGGTGGACTACG
 GCTTCGGGCTGAGCCTGGAGGCTGCACGCCACCATGAGAACAGCTCTCGGGGCGCAGTGGTGGGCGTGGG
 CGACGAGGGCCGCGCCCTGCCGATGTCATGCAGGGCATGGTGTCTCAGCTCCATGCAGCACTTCAGCGAG
 GCCTTCCACCAGGTCTGGGCGAAAAGCATAAGCGTAGCCACCTGGCCGAGGCCGAGGGCCATAGGGACA
 CTTGTGATGAAGACTCGGTGGCCGGTGAAGTACAGCCGATAGACGATGGCACTGTTAATGGTCGTGGCTG
 CTCCCCCGGCGAATCGGCTTCGGGGGCTGTGCCAAAAGCTGCTGCTGGGTAGCCCCAGCTCGCTGAGC
 CCCTTCTCCAAGCGCATCAAGCTGGAGAAGGAGTTTACCTGCCCCGGCCGCGATGCCTAACACGGAGA
 ACGTGTATTGCGAGTGGCTCGCTGGCTATGCGGCCTCCAGGCAGCTCAAAGATCCCTTCTTACTTTCGG
 AGACTCCAGACAAATCGCCTTTTGCTCTCATCAGAGCACTCTCGGAGAACGGGAGCTTGCCTTCTCC
 ACACCGCCCGGGAGCTGGACGGAGGGATCTCAGGGCGCAGCGGCACAGGAAGTGGAGGGAGCACGCCCC
 ATATTAGTGGTCCGGGCCCGGCAGGCCAGCTCAAAGAGGGCAGACGACGACACTTGTCTTCACA
 CACCCCGTTTCGGCGTAGTACCCCGCAGCTCAAGATGTGTGGCAGTTTTTCGGATGGAAGCTCAAGAACC
 CTTAAGTTC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR221882 representing NM_016707
 Red=Cloning site Green=Tags(s)


MSRRKQGKPKHL SKREFSPEPLEAILTDDEPDHGPLGAPEGDHDL TCGQCQMNFLGDILIFIEHKRKQ
 CNGSLCLEKGVDPKPPSPSPIEMKKASNPVEVGIVTPEDDDCLSTSSRGICPKQEHADKLLHWRGLSSP
 RSAHGALIPTPGMSAEYAPQGICKDEPSSYTCTTCKQPFTSAWFLQHAQNTHGLRIYLESEHGSP LTPR
 VGIPSGLGAECPSQPPLHGIHIADNNPFNLLRIPGSVSREASGLAEGRFPTPLFSPPPRHLDPHRIE
 RLGAEEMLATHHPSAFDRVLRNLNPMAMEPPAMDFSRRLRELAGNTSSPPLSPGRPSMQRLLPFQPGS
 KPPFLATPPLPPLQSAPPPSQPPVKSKEFCGKTFKFQSNLVVHRRSHTGEKPYKCNLCDHACTQASKL
 KRHMKTHMHKSSPMTVKSDGLSTASSPEPGTSDLVGSASSALKSVVAKFKSENDPNLIPENGDEEEEEE
 DEEEEEEEEEEEELTESERVDYGFGLSLEAARHHENSSRGAVVGVDGRALPDVMQGMVLSMQHFSE
 AFHQVLGEKHKRSHLAEAEGRDTCDEDSVAGESDRIDDGTVNGRGCSPGESASGGLSKKLLLGSPSSLS
 PFSKRIKLEKEFDLPPAAMPNTENVYSQWLAGYAASRQLKDPFLTFGDSRQSPFASSEHSSENGSLRFS
 TPPGELDGGISGRSGTGGGSTPHISGPGPRPSSKEGRRSDTCPSTPVRRTSTPRAQDVWQFSDGSSRT
 LKF

SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9094_d08.zip

Restriction Sites: Sgfl-RsrII

Cloning Scheme:



ACCN: NM_016707


ORF Size: 2319 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_016707.3, NP_057916.1</u>
RefSeq Size:	3425 bp
RefSeq ORF:	2322 bp
Locus ID:	14025
UniProt ID:	<u>Q9QYE3</u>
Cytogenetics:	11 A3.2
MW:	83.9 kDa
Gene Summary:	Functions as a myeloid and B-cell proto-oncogene. May play important roles in leukemogenesis and hematopoiesis. An essential factor in lymphopoiesis, is required for B-cell formation in fetal liver. May function as a modulator of the transcriptional repression activity of ARP1.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR221882