

Product datasheet for **MR221888A1V**

Mouse Bcl11a (NM_001159290) AAV Particle

Product data:

Product Type:	AAV Particles
Tag:	Myc-DDK
Symbol:	Bcl11a
Synonyms:	2810047E18Rik; BCL-11A; Ctbp1; D930021L15Rik; Evi9; Evi9a; Evi9b; Evi9c; mKIAA1809
Mammalian Cell	None
Selection:	
Vector:	pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >MR221888 representing NM_001159290
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGAATTTCCCACTGGGGACATTCTTATTTTATCGAGCACAAACGGAACAATGCAATGGCAGCCTCT
GCTTAGAAAAAGGTGGATAAGCCGCTTCCCTTCTCCATCGAGATGAAAAAGGCATCCAATCCTGT
GGAGGTTGGCATCCAGGTCACGCCAGAGGATGACGATTGTTTATCAACGTCATCTAGAGGAATTTGCCCC
AAACAGGAACATAGCAGATAAACTTCTGCACTGGAGGGGCTGTCTCTCCTCGGTCTGCACACGGAG
CTCTAATCCCAACGCGGGATGAGTGCAGAAATATGCCCGCAGGGTATTTGTAAGATGAGCCAGCAG
CTACACATGTACAACCTTGCAACAGCCATTACCAAGTGCATGGTTTCTCTTGAACACGCACAGAACACT
CATGGATTAAGAACTACTTAGAAAGTGAACACGGAAGTCCCTGACCCGCGGGTCTTCACACACCCC
CGTTCGGCGTAGTACCCGCGAGCTCAAGATGTGTGGCAGTTTTCGGATGGAAGCTCAAGAACCCTTAAG
TTCTGAGAACTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGAGATAAGGTTTAA

Protein Sequence: >MR221888 representing NM_001159290
Red=Cloning site Green=Tags(s)

MNFPGLDILIFIEHKRKQCNGSLCLEKGVDPKPPSPPIEMKKASNPVEVGIVQTPEDDDCLSTSSRGICP
KQEHIAADKLLHWRGLSSPRSAHGALIPTPGMSAEYAPQGICKDEPSSYTCTCKQPFSAWFLQHAQNT
HGLRIYLESEHGSPLTPRVLHTPPFGVVPRELKMGCSFRMEAQPLSSEKL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse



Serotype:	AAV-2
ACCN:	NM_001159290
ORF Size:	573 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001159290.1 , NP_001152762.1
RefSeq Size:	1480 bp
RefSeq ORF:	576 bp
Locus ID:	14025
UniProt ID:	Q9QYE3
Cytogenetics:	11 A3.2
MW:	21.5 kDa