

Product datasheet for **MR221898L3V**

Pfkm (NM_001163488) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pfkm (NM_001163488) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pfkm
Synonyms:	A1131669; ATP-PFK; Pfk-4; PFK-A; PFK-M; Pfk4; PfkA; Pfkx
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001163488
ORF Size:	2343 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221898).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001163488.1 , NP_001156960.1
RefSeq Size:	2830 bp
RefSeq ORF:	2343 bp
Locus ID:	18642
UniProt ID:	P47857
Cytogenetics:	15 F1
Gene Summary:	Catalyzes the phosphorylation of D-fructose 6-phosphate to fructose 1,6-bisphosphate by ATP, the first committing step of glycolysis.[UniProtKB/Swiss-Prot Function]



[View online »](#)