

## Product datasheet for **MR221909L4V**

### Tnik (NM\_001163008) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tnik (NM_001163008) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tnik
Synonyms:	1500031A17Rik; 4831440119Rik; AI451411; C530008O15Rik; C630040K21Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001163008
ORF Size:	3969 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221909).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001163008.1</a> , <a href="#">NP_001156480.1</a>
RefSeq Size:	7035 bp
RefSeq ORF:	3972 bp
Locus ID:	665113
UniProt ID:	<a href="#">P83510</a>
Cytogenetics:	3 A3



[View online »](#)

**Gene Summary:**

Serine/threonine kinase that acts as an essential activator of the Wnt signaling pathway. Recruited to promoters of Wnt target genes and required to activate their expression. May act by phosphorylating TCF4/TCF7L2. Appears to act upstream of the JUN N-terminal pathway. May play a role in the response to environmental stress. Part of a signaling complex composed of NEDD4, RAP2A and TNIK which regulates neuronal dendrite extension and arborization during development. More generally, it may play a role in cytoskeletal rearrangements and regulate cell spreading (By similarity).[UniProtKB/Swiss-Prot Function]