

Product datasheet for MR221926

Syt7 (NM_018801) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Syt7 (NM_018801) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Syt7
Synonyms: AI851541; B230112P13Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR221926 representing NM_018801
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTACCGGGACCCGGAGGCGGCCAGCCAGGGGCACCTACCGCGATGTCCTGCTGGTCTCTGCAATCA
TCACCGTCAGCCTTAGCGTCACTATCGTCTCTGCGGCCGTGCCACTGGTGTGACGCGAACTGGGCAA
ACGCTACAAGAATTCCTTGGAGACGGTGGGCACGCCAGACTCGGGCGTGGGCGCGGTGAGAAGAAAGCC
ATCAAGTTGCCTGCAGGAGGGAAGGCTGTGAATACAGCCCAAGTCCCGGCCAGACGCCACACGATGAGT
CTGACCGCAGAACGGAGACCCGTTCTCTGTCTCGGACCTCGTCAACTCCCTTACCAGCGAGATGCTCAT
GCTCTCCCGGGTTCGAGGAGGATGAGGCCACAGGGGCTGCAGCCGAGAGAACCTGGGCCGAATCCAG
TTCAGTGTGGCTACAACCTTCCAAGAGTCCACACTCACCGTGAAGGTCATGAAGGCCAAGAGCTGCCAG
CCAAGGACTTCAGCGGTACTAGTGACCCCTTTGTCAAGATCTACCTGCTACCCGACAAGAAGCACAACT
GGAGACCAAGGTGAAGCGGAAGAATCTAAACCCGACTGGAATGAGACCTTTCTATTTGAAGGGTTTCCC
TACGAGAAAGTGGTGCAGAGGGTCTCTACCTCCAGGTCTGGATTATGACCGTTTCAGCCGCAATGACC
CCATTGGGAGGTGTCCATCCCTCTGAACAAGTGGACCTGACCCAGATGCAGACCTTCTGGAAGGATCT
GAAGCCATGCAGCGATGGGAGTGGGAGCCGAGGGGAGCTGCTCTGTCCCTCTGCTACAACCCCTCTGCC
AACTCCATCATCGTGAACATCATCAAAGCTCGAAACCTCAAAGCCATGGACATCGGGGGCACATCAGACC
CCTATGTGAAGGTGTGGCTGATGTATAAAGACAAGCGGGTAGAGAAAAAGAACCGTGACAAGAAGAG
GAACCTGAACCCATCTTCAATGAGTCTTTCGCTTCGACATACCCACGAGAAGCTGAGGGAGACCACG
ATCATCATCACTGTCATGGACAAAGACAAGCTCAGCCGCAACGACGTCATCGGCAAGATCTACCTGTCT
GGAAGAGCGGACCAGGTGAAGTGAACACTGGAAGGACATGATCGCTCGTCCCGGCAGCCTGTGGCCCA
GTGGCACCAGCTGAAAGCC

ACGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online >](#)

Protein Sequence: >MR221926 representing NM_018801
Red=Cloning site Green=Tags(s)

MYRDPEAASPGAPTRDVLLVSAIITVSLSVTIVLCGLCHWCQRKLGKRYKNSLETVGTDPDSGRGRGEKKA
 IKLPAGGKAVNTAPVPGQTPHDESDRRTETRSSVSDLVNSLTSEMLMLSPGSEEDAHEGCSRENLGRIQ
 FSVGYNFQESTLTVKVMKAQELPAKDFSGTSDPFVKIYLLPDKKHKLETKVKKRNLNPHWNETFLFEGFP
 YEKVVQRVLYLQVLDYDRFSRNDPIGEVSIPLNKVDLTQMOTFWKDLKPCSDGSGSRGELLLSLCYNPSA
 NSIIVNIIKARNLKAMDIGGTSDPYVKVWLMYKDKRVEKKTIVTKRNLNPIFNESFAFDIPTEKLRRETT
 IITVMDKDKLSRNDVIGKIYLSWKS GPGEVKHWKDMIARPRQPV AQWHQLKA


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9007_c02.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_018801

ORF Size: 1209 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_018801.3](#), [NP_061271.1](#)

RefSeq Size: 6342 bp

RefSeq ORF: 1212 bp

Locus ID: 54525


UniProt ID: [Q9R0N7](#)

Cytogenetics: 19 A

MW: 45.9 kDa

Gene Summary:

Ca(2+) sensor involved in Ca(2+)-dependent exocytosis of secretory and synaptic vesicles through Ca(2+) and phospholipid binding to the C2 domain. Ca(2+) induces binding of the C2-domains to phospholipid membranes and to assembled SNARE-complexes; both actions contribute to triggering exocytosis. SYT7 binds Ca(2+) with high affinity and slow kinetics compared to other synaptotagmins (PubMed:26738595). Involved in Ca(2+)-triggered lysosomal exocytosis, a major component of the plasma membrane repair (By similarity). Ca(2+)-regulated delivery of lysosomal membranes to the cell surface is also involved in the phagocytic uptake of particles by macrophages (PubMed:16982801, PubMed:21041449). Ca(2+)-triggered lysosomal exocytosis also plays a role in bone remodeling by regulating secretory pathways in osteoclasts and osteoblasts (PubMed:18539119). Involved in cholesterol transport from lysosome to peroxisome by promoting membrane contacts between lysosomes and peroxisomes: probably acts by promoting vesicle fusion by binding phosphatidylinositol-4,5-bisphosphate on peroxisomal membranes (PubMed:25860611). Acts as a key mediator of synaptic facilitation, a process also named short-term synaptic potentiation: synaptic facilitation takes place at synapses with a low initial release probability and is caused by influx of Ca(2+) into the axon terminal after spike generation, increasing the release probability of neurotransmitters (PubMed:24569478, PubMed:26738595). Probably mediates synaptic facilitation by directly increasing the probability of release (PubMed:26738595). May also contribute to synaptic facilitation by regulating synaptic vesicle replenishment, a process required to ensure that synaptic vesicles are ready for the arrival of the next action potential: SYT7 is required for synaptic vesicle replenishment by acting as a sensor for Ca(2+) and by forming a complex with calmodulin (PubMed:24569478). Also acts as a regulator of Ca(2+)-dependent insulin and glucagon secretion in beta-cells (PubMed:18308938, PubMed:19171650). Triggers exocytosis by promoting fusion pore opening and fusion pore expansion in chromaffin cells (PubMed:20956309). Also regulates the secretion of some non-synaptic secretory granules of specialized cells (By similarity). [UniProtKB/Swiss-Prot Function]

Product images:


Circular map for MR221926