

## Product datasheet for **MR221942L4V**

### MgII (NM\_001166250) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	MgII (NM_001166250) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	MgII
Synonyms:	AA589436; Ma; MagI; Mgl
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001166250
ORF Size:	396 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221942).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001166250.1</a> , <a href="#">NP_001159722.1</a>
RefSeq Size:	1204 bp
RefSeq ORF:	399 bp
Locus ID:	23945
UniProt ID:	<a href="#">O35678</a>
Cytogenetics:	6 D1



[View online »](#)

**Gene Summary:**

This gene encodes a monoglyceride lipase, which catalyzes the hydrolysis of monoglycerides into fatty acids and glycerol. This enzyme is also thought to hydrolyze the endocannabinoid 2-arachidonoylglycerol. Alternatively spliced transcript variants have been described. [provided by RefSeq, Oct 2009]