

## Product datasheet for **MR221946**

### Sel1l (NM\_001039089) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Sel1l (NM_001039089) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Sel1l
Synonyms:	AW493766; mKIAA4137; Sel1h
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>MR221946 representing NM\_001039089  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**

ATGCGAGTCCGCGTCAGGCTGTCTGCTGCTCTGCGCGGTGCTCCTGGGCTCGGCAGCCGCGACCT  
 CGGATGACAAAACCTAACCCAGGATGACTCCTTAGATTCCAAGAGTTCTTTGCCACAGATGAGTCAGTGAA  
 GGACCACACCACCACGGGCAAAGTAGTTGCTGGCCAGATATTTGTTGATTCTGAAGAAGCAGAAAGTGAA  
 TCCCTTCTTCAGGACGAGGAAGTAGCTCCAAGACCCAGGAGGAAGAGATCAGCTTTTTAGAATCTCCGA  
 ATCCAAGCAGCAAGACCTACGAAGAATAAAGAGAGTGCAGGAGCCAGTCTTGACTGCCATTGAAGGTAC  
 GGCGCACGGGAGCCCTGCCACTTCCCTTTCTTTTCTGGATAAGGAGTATGATGAGTGCACCTCAGAC  
 GGGAGGGAAGATGGCAGACTGTGGTGTGCCACAACCTATGACTACAAGACAGATGAGAAGTGGGGCTTCT  
 GCGAAACTGAAGAAGATGCTGCCAAAAGACGACAGATGCAGGAAGCAGAGATGATCTATCAGGCCGGAT  
 GAAGATACTGAATGGAAGCAATAGGAAGAGCCAAAAGAGAGAAGCATATCGGTACCTTCAGAAGGCAGCA  
 GGCATGAATCACACAAAAGCCCTGGAGAGAGTGTCTATGCTCTCTTGTGGTATTACCTCACACAGA  
 ATATCCAGGCAGCCAAAGAGATGTTGAGAAAAGTACTGAGGAAGGGTCTCCCAAAGGACAGACTGGTCT  
 TGGCTTCTCTACGCTTCTGGGCTTGGTGTAAATCAAGTCAGGCAAAGGCTCTTGATATATTACTTTT  
 GGAGCTCTGGAGGCAACCTGATAGCCCATATGATTTTGGGTTACCGCTACTGGGCTGGCATCGGAGTCC  
 TCCAGAGTTGTGAGTCCGCACTGACCCATTATCGTCTTGTGGCAATCATGTTGCTAGTGATATCTCCCT  
 AACTGGAGGCTCTGTAGTCCAGAGAATACGGCTGCCCGATGAAGTGGAAAACCCGGGGATGAACAGTGGG  
 ATGCTGGAAGAAGACCTGATTCACTATTACCAGTTCCTAGCTGAGAAGGGTGACGTCCAAGCACAGTTG  
 GCCTGGGACAGCTGCATCTGCATGGAGGGCGTGGAGTAGAACAAGAAATCACCAGAGAGCGTTTACTT  
 CAACTTAGCAGCAAATGCTGGCAATTACATGCTATGGCCTTCTGGGAAAGATGATTCTGAAGGAAGT  
 GACATCGTACCTCAGAGTAAAGAGACGGCACTTCACTACTTTAAGAAAGCTGCTGACATGGGCAACCCCG  
 TGGGACAGAGCGGGCTTGGAAATGGCCTACCTCTACGGAAGAGGCGTTCAAGTTAATTATGACCTGGCCCT  
 CAAGTATTTCCAGAAAGCTGCTGAGCAAGGCTGGGTGGACGGGACAGTGCAGCTCGGCTCTATGTACTAC  
 AATGGCATTGGAGTCAAGAGAGATTATAAGCAGGCCTTGAAGTATTTAATCTGGCTTCTCAAGGAGGCC  
 ATATCTTGGCTTCTATAACCTCGCACAGATGCACGCCAGCGGCACAGGGGTGATGCGGTCTGTACAC  
 TGCAGTGGAGTTGTTAAGAATGTGTGAGCGAGGTGCTGGTCCAGAGAGACTGATGACTGCCTACAAC  
 AGCTATAAGGATGAGGACTACAATGCTGCAGTGGTCCAGTACCTCCTGCTGGCTGAGCAGGGCTACGAGG  
 TGGCGCAGAGCAACGCAGCCTTATCCTCGACCAGAGAGAAGCAACCATTGTAGGTGAGAATGAACTTA  
 CCCCAGAGCTTTACTGCATTGGAACAGGGCCGCTCCCAAGGTTACACTGTGGCTAGAATTAAGCTTGG  
 GACTACCCTTCTATGGCTTTGGCACTGATGTGGATTATGAGACCGCATTATTACCTACCGCTGGCTT  
 CTGAGCAGCAGCAGCGCCCAAGCTATGTTTAACTGGGCTACATGCACGAGAAGGGCTAGGCATTAA  
 ACAGGACATTACCTTGCAAAACGCTTTTATGACATGGCAGCCGAAGCTAGCCAGATGCACAAGTACCT  
 GTGTTCTCGCACTCTGCAAAATAGGTGTCGTCTATTTCTTACAGTACATACGGGAAGCAAAACATTCGAG  
 ATCTATTACACAACCTGGATATGGACCAGCTTTTGGGACCCGAGTGGGACCTTTACCTCATGACCATCAT  
 TGCAGTCTCTTGGGTACAGTCATAGCTTACAGGCAGCGACAGCACCAGGACATACCAGTTCACAGGCC  
 CCAGGGCCACGGCCGGCTCCTCCCAGCAGGAAGGACCACCAGAGCAGCAGCCACCACAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR221946 representing NM\_001039089  
Red=Cloning site Green=Tags(s)


MQVRVRLSLLLLCAVLLGSAAATSDDKTNQDDSLDSKSSLPTDESVKDHTTTGKVVAGQIFVDSEEA EVE  
SLLQDEEDSSKTQEEIISFLESPNPSSKTYEELKRVKPVLTAIEGTAHGEPCHFPLFLDKEYDECTSD  
GREDGRLWCATTYDYKTDEKWGFCETEEDAARRQMQEAEMIYQAGMKILNGSNRKSQKREAYRYLQKAA  
GMNHTKALERSYALLFGDYLTQNIQAAKEMFEKLTEEGSPKGGTGLGFLYASGLGVNSSQAKALVYYTF  
GALGGNLI AHMILGYRYWAGIGVLQSCESALTHYRLVANHVASDISLTGGSVVQRIRLPDEVENPGMNSG  
MLEEDLIQYYQFLAEKGDVQAQVGLGQLHLHGGRGVEQNHQRAFDYFNLAANAGNSHAMAFLGKMYSEGS  
DIVPQSNETALHYFKKAADMGNPVGQSGLGMAYL YGRGVQVNYDLALKYFQKAAEQGWVDGQLQLGSMYY  
NGIGVKRDYKQALKYFNLASQGGHILAFYNLAQMHASGTGVMRSCHTAVELFKNVCERGRWSERLMTAYN  
SYKDEDYNAAVVQYLLAEQGYEVAQSNAAFILDQREATIVGENETYPRALLHWNRAASQGYTVARIKLG  
DYHFYGFGTDDVDYETAFIHYRLASEQQHSAQAMFNLGYMHEKGLGIKQDIHLAKRFYDMAAEASPDQVP  
VFLALCKLGVVYFLQYIREANIRDLFTQLDMDQLLGPWDLYLMTIALLLGTVIAYRQRQHQDIPVPRP  
PGPRPAPPQQEGPPEQQPPQ

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9004\\_a10.zip](https://cdn.origene.com/chromatograms/mm9004_a10.zip)

**Restriction Sites:** Sgfl-Mlul

Cloning Scheme:



ACCN: NM\_001039089

ORF Size: 2370 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001039089.1](#), [NP\\_001034178.1](#)

**RefSeq Size:** 6300 bp

**RefSeq ORF:** 2373 bp

**Locus ID:** 20338


**UniProt ID:** [Q9Z2G6](#)

**Cytogenetics:** 12 D3

**MW:** 88.8 kDa

**Gene Summary:** Plays a role in the endoplasmic reticulum quality control (ERQC) system also called ER-associated degradation (ERAD) involved in ubiquitin-dependent degradation of misfolded endoplasmic reticulum proteins (PubMed:25066055, PubMed:24453213). Enhances SYVN1 stability (PubMed:24453213). Plays a role in LPL maturation and secretion (PubMed:25066055). Required for normal differentiation of the pancreas epithelium, and for normal exocrine function and survival of pancreatic cells (PubMed:20170518, PubMed:24453213). May play a role in Notch signaling (PubMed:20170518).[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR221946