

## Product datasheet for **MR221951L3V**

### Abcc6 (NM\_018795) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Abcc6 (NM_018795) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Abcc6
Synonyms:	Abcc1b; DCC; Dysca; dyscalc; Dyscalc1; Mr; Mrp6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_018795
ORF Size:	4494 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221951).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_018795.2</a> , <a href="#">NP_061265.2</a>
RefSeq Size:	4958 bp
RefSeq ORF:	4497 bp
Locus ID:	27421
UniProt ID:	<a href="#">Q9R1S7</a>
Cytogenetics:	7 B3



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the MRP subfamily which is involved in multi-drug resistance. The specific function of this protein is unknown; however, a similar rat protein has been identified as the major canalicular bile salt export pump of liver. [provided by RefSeq, Jul 2008]