

Product datasheet for **MR221972L3V**

Fgf12 (NM_183064) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fgf12 (NM_183064) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fgf12
Synonyms:	AV114868; B230343J05Rik; FGF-12; FHF-1; Fhf1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_183064
ORF Size:	732 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221972).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_183064.3 , NP_898887.1
RefSeq Size:	5449 bp
RefSeq ORF:	732 bp
Locus ID:	14167
UniProt ID:	P61329
Cytogenetics:	16 19.38 cM
Gene Summary:	Involved in nervous system development and function. Promote neuronal excitability by elevating the voltage dependence of neuronal sodium channel SCN8A fast inactivation. [UniProtKB/Swiss-Prot Function]



[View online »](#)