

Product datasheet for MR221993

Abhd13 (NM_001081119) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Abhd13 (NM_001081119) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Abhd13
Synonyms:	1110065L07Rik; AI463703; AI788994
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR221993 representing NM_001081119 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAAAAATCCTGGATGCTGTGGAGCTTCATTGAAAGATGGCTACTGGCCTTGGCGTCTGGTCTTGGG
CTCTCTGCCGCATCTCTGTACCTTTGATAGTGACTTTTCATCTGTATGGAGGCATTGTCTACTTTT
GTTAATATTCGTATCAATTGCAGGTATCTATATAAATCCAAGATGCTTGCTCTACTTTCCAGAACAG
CCATCTTCTTCCCGCCTTATGTTCCCATGCCACGGGTATCCCATGAAAACATTTTCATCAGAACCA
AAGATGGAGTGGCCTGAATCTGATCTTGGTGAGATACTGGAGACAATCCCCCTACTGCCGACCAT
AATTTATTTCCATGGAAATGCCGCAACATAGGCCACAGGTTACCAAACGCACTGCTAATGCTGGTGAAC
CTCAGAGTGAATCTTGTGCTTGTGGATTACCGGGGTACGGGAAGAGCGAAGGAGAGGCCAGTGAAGAAG
GGCTGTATTTGGATTCTGAAGCTGTGCTAGACTATGTCATGACGAGACCCGACCTCGATAAGACCAAGGT
TTTCTTTTGGCCGCTCCCTGGGAGGAGCAGTGGCCATTCACCTTGGCTTCTGAGAATTCTCACAGGATT
TCAGCCATCATGGTGGAGAACACATTCCTAAGCATACCACATGGCCAGCACTATTTTCATTCTTCC
CAATGCGTTACCTTCTTATGGTGCTATAAAAATAAGTTTCTGTCTACAGAAAGATCTCCAGTGCAG
GATGCCTTCTCTTTCATCTCTGGCCTCTTGACCAGTTGATCCACCCGTGATGATGAAACAGCTCTAC
GAGCTCTCCATCCCGGACTAAGAGATTAGCCATTTTCCAGATGGGACGCACAACGACACATGGCAGT
GCCAAGGCTACTTCACTGCACTAGAACAATTTATCAAGGAAGTGATCAAGAGTCATTCTCCGAAGACAT
GACAAAAACATCATCCAATGTCACCATCATC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR221993 representing NM_001081119
Red=Cloning site Green=Tags(s)

MEKSWMLWSFIERWLLALASWSWALCRISLLPLIVTFHLYGGIVLLLLIFVSIAGILYKFQDVLLYFPEQ
 PSSSRLVYPMPTGIPHENIFIRTKDGVRLNLILVRYTGDNSPYCPTIIYFHGNAGNIGHRLPNALLMLVN
 LRVNLVLVDYRGYKSEGEASEEGLYLDSEAVLDYVMTRPDLDKTKVFLFGRSLGGAVAIHLASENSHRI
 SAIMVENTFLSIPHMASTLFSFFPMRYLPLWCYKNKFLSYRKISQCRMPSLFISGLSDQLIPPVMMKQLY
 ELSPSRTKRLAIFPDGTHNDTWQCQGYFTALEQFIKEVIKSHSPEDMTKTSNVITII

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001081119

ORF Size: 1011 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001081119.1](#), [NP_001074588.1](#)

RefSeq Size: 5074 bp

RefSeq ORF: 1014 bp


Locus ID: 68904

UniProt ID: [Q80UX8](#)

Cytogenetics: 8 A1.1

MW: 38.5 kDa

Product images:



Circular map for MR221993