

Product datasheet for **MR222010L3V**

Palm3 (NM_028877) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Palm3 (NM_028877) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Palm3
Synonyms:	4432412L15Rik; AL024300
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_028877
ORF Size:	2202 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222010).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_028877.1 , NP_083153.1
RefSeq Size:	2449 bp
RefSeq ORF:	2205 bp
Locus ID:	74337
UniProt ID:	A2TIV2
Cytogenetics:	8 C2
Gene Summary:	ATP-binding protein, which may act as a adapter in the Toll-like receptor (TLR) signaling. [UniProtKB/Swiss-Prot Function]



[View online »](#)