

Product datasheet for **MR222022L3V**

Fetub (NM_001083904) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fetub (NM_001083904) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fetub
Synonyms:	2310011O17Rik; AI255764; D17980; Gugu
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001083904
ORF Size:	924 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222022).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001083904.1 , NP_001077373.1
RefSeq Size:	1365 bp
RefSeq ORF:	927 bp
Locus ID:	59083
UniProt ID:	Q9QXC1
Cytogenetics:	16 13.76 cM



[View online »](#)

Gene Summary:

Protease inhibitor required for egg fertilization. Required to prevent premature zona pellucida hardening before fertilization, probably by inhibiting the protease activity of ASTL, a protease that mediates the cleavage of ZP2 and triggers zona pellucida hardening. [UniProtKB/Swiss-Prot Function]