

Product datasheet for **MR222028L4V**

Tmem167 (NM_025335) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Tmem167
Synonyms:	0610041E09Rik; 5730424F14Rik; AU041184; AW537806; ENSMUSG00000060669; Gm10085; Tmem167a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025335
ORF Size:	216 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222028).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025335.3 , NP_079611.2
RefSeq Size:	10809 bp
RefSeq ORF:	219 bp
Locus ID:	66074
UniProt ID:	Q9CR64
Cytogenetics:	13 C3



Gene Summary:

Involved in the early part of the secretory pathway.[UniProtKB/Swiss-Prot Function]