

## Product datasheet for MR222053

### Ap3s1 (NM\_009681) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ap3s1 (NM\_009681) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Ap3s1  
**Synonyms:** [s]  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR222053 representing NM\_009681  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGATCAAGGCCATCCTCATCTTCAACAACCACGGGAAGCCGGCTCTCCAAGTTCTACCAGCCCTATA  
 GTGAAGACACGCAACAGCAAATCATCAGGGAGACTTCCATTTGGTGTCTAAGCGCGATGAGAACGTTTG  
 TAATTTCTAGAAGGAGGATTATTAATTGGAGGCTCTGACAACAAGCTCATTTACAGACATTATGCAACA  
 CTATATTTTGTCTTCTGTGTGGACTCCTCAGAAAGTGAAGTGGCATTATAGATCTAATTCAGTATTTG  
 TGGAAACATTAGACAAATGTTTTGAAAATGTTTGTGAAGTGGATTAATATTCCATGTAGACAAGGTTCA  
 TAATATTCTTGCAGAAATGGTGTATGGGGGAATGGTATTGGAGACCAACATGAATGAGATTGTACACAA  
 ATTGATGCACAAAATAAACTGGAGAAATCTGAGGCTGGCTTAGCAGGAGCTCCAGCCCGTGTGTATCAG  
 CTGTAAAGAATATGAATCTTCTGAGATCCCAAGAAATATTAACATTGGTACATCAGTATAAAAGTGCC  
 AAACCTGCCCTCTTTAAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR222053 representing NM\_009681  
 Red=Cloning site Green=Tags(s)

MIKAILIFNNHGKPRLSKFYQPYSEDYQQIIRETFHLVSKRDENVCFLEGGLLIGGSDNKLIRHYAT  
 LYFVFCVDSSESELGILDLIQVFVETLDKCFENVCELDLIFHVDKVHNI LAEMVMGGMVLETNMNEIVTQ  
 IDAQNKLEKSEAGLAGAPARAVSAVKNMNLP EIPRNINIGDISIKVPLPSFK

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**



**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_009681

**ORF Size:** 579 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_009681.5](#), [NP\\_033811.1](#)

**RefSeq Size:** 1479 bp

**RefSeq ORF:** 582 bp

**Locus ID:** 11777


**UniProt ID:** [Q9DCR2](#)

**Cytogenetics:** 18 C

**MW:** 22.2 kDa

**Gene Summary:** This gene encodes the sigma subunit of the heterotetrameric adaptor protein complex AP-3 which is involved in the formation of specialized lysosome-related compartments such as melanosomes. Alternate splicing of this gene results in multiple transcript variants encoding different isoforms. Pseudogenes of this gene are found on chromosomes 1, 8, 16, 17 and X. [provided by RefSeq, Dec 2014]

**Product images:**



Circular map for MR222053