

Product datasheet for **MR222073L4V**

Fga (NM_001111048) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fga (NM_001111048) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fga
Synonyms:	Fi; Fib
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001111048
ORF Size:	2367 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222073).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001111048.1 , NP_001104518.1
RefSeq Size:	2652 bp
RefSeq ORF:	2370 bp
Locus ID:	14161
UniProt ID:	E9PV24
Cytogenetics:	3 36.96 cM



[View online »](#)

Gene Summary:

This gene encodes a subunit of the coagulation factor fibrinogen, which is a component of the blood clot. The encoded protein is proteolytically processed by thrombin during the conversion of fibrinogen to fibrin. Mice lacking the encoded protein display bleeding in the peritoneal cavity, skin and soft tissues around joints immediately after birth, and are predisposed to spontaneous fatal abdominal hemorrhage as they grow. Pregnant mice lacking the encoded protein succumb to uterine bleeding during gestation. Alternative splicing results in multiple transcript variants encoding different isoforms that may undergo similar proteolytic processing. [provided by RefSeq, Nov 2015]