

## Product datasheet for **MR222093L3V**

### **Pitx3 (NM\_008852) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

|                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | Lentiviral Particles                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Product Name:             | Pitx3 (NM_008852) Mouse Tagged ORF Clone Lentiviral Particle                                                                                                                                                                                                                                                                                                                                                                                       |
| Symbol:                   | Pitx3                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Synonyms:                 | ak; Ptx3                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Mammalian Cell Selection: | Puromycin                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)                                                                                                                                                                                                                                                                                                                                                                                                               |
| Tag:                      | Myc-DDK                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| ACCN:                     | NM_008852                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| ORF Size:                 | 906 bp                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(MR222093).                                                                                                                                                                                                                                                                                                                                                                                     |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.                                                                                                                                                                                                                                                                                                               |
| RefSeq:                   | <a href="#">NM_008852.4</a> , <a href="#">NP_032878.1</a>                                                                                                                                                                                                                                                                                                                                                                                          |
| RefSeq Size:              | 1379 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| RefSeq ORF:               | 909 bp                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Locus ID:                 | 18742                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| UniProt ID:               | <a href="#">O35160</a>                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Cytogenetics:             | 19 38.75 cM                                                                                                                                                                                                                                                                                                                                                                                                                                        |


[View online »](#)

**Gene Summary:**

Transcriptional regulator which is important for the differentiation and maintenance of meso-diencephalic dopaminergic (mdDA) neurons during development. In addition to its importance during development, it also has roles in the long-term survival and maintenance of the mdDA neurons. Activates NR4A2/NURR1-mediated transcription of genes such as SLC6A3, SLC18A2, TH and DRD2 which are essential for development of mdDA neurons. Acts by decreasing the interaction of NR4A2/NURR1 with the corepressor NCOR2/SMRT which acts through histone deacetylases (HDACs) to keep promoters of NR4A2/NURR1 target genes in a repressed deacetylated state. Essential for the normal lens development and differentiation. Plays a critical role in the maintenance of mitotic activity of lens epithelial cells, fiber cell differentiation and in the control of the temporal and spatial activation of fiber cell-specific crystallins. Positively regulates FOXE3 expression and negatively regulates PROX1 in the anterior lens epithelium, preventing activation of CDKN1B/P27Kip1 and CDKN1C/P57Kip2 and thus maintains lens epithelial cells in cell cycle.[UniProtKB/Swiss-Prot Function]