

## Product datasheet for **MR222101L3V**

### Adgrg2 (NM\_178712) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Adgrg2 (NM_178712) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Adgrg2
Synonyms:	AW212196; B830041D06Rik; Gpr64; Me6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_178712
ORF Size:	3027 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222101).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_178712.4</a> , <a href="#">NP_848827.1</a>
RefSeq Size:	4727 bp
RefSeq ORF:	3030 bp
Locus ID:	237175
UniProt ID:	<a href="#">Q8CJ12</a>
Cytogenetics:	X F4
Gene Summary:	Orphan receptor. Could be involved in a signal transduction pathway controlling epididymal function and male fertility. May regulate fluid exchange within epididymis.[UniProtKB/Swiss-Prot Function]



[View online »](#)