

## Product datasheet for **MR222112L3V**

### **Bex2 (NM\_009749) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Bex2 (NM_009749) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Bex2
Synonyms:	AL024066; Bex1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_009749
ORF Size:	387 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222112).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_009749.2</a> , <a href="#">NP_033879.1</a>
RefSeq Size:	902 bp
RefSeq ORF:	390 bp
Locus ID:	12069
UniProt ID:	<a href="#">Q9WTZ8</a>
Cytogenetics:	X 57.37 cM
Gene Summary:	Regulator of mitochondrial apoptosis and G1 cell cycle. Regulates the level of PP2A regulatory subunit B and PP2A phosphatase activity (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)