

Product datasheet for MR222123L4

Ghr (NM_001048147) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ghr (NM_001048147) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Ghr
Synonyms:	AA986417; GHBP; GHR/BP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222123).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

Locus ID: 14600

Cytogenetics: 15 1.84 cM

Gene Summary: Receptor for pituitary gland growth hormone involved in regulating postnatal body growth. On ligand binding, couples to, and activates the JAK2/STAT5 pathway (By similarity).
[UniProtKB/Swiss-Prot Function]