

## Product datasheet for **MR222171L3V**

### Actn2 (NM\_033268) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Actn2 (NM_033268) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Actn2
Synonyms:	1110008F24Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_033268
ORF Size:	2685 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222171).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_033268.4</a>
RefSeq Size:	3091 bp
RefSeq ORF:	2685 bp
Locus ID:	11472
UniProt ID:	<a href="#">Q9J191</a>
Cytogenetics:	13 4.54 cM
Gene Summary:	F-actin cross-linking protein which is thought to anchor actin to a variety of intracellular structures. This is a bundling protein (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)