

Product datasheet for MR222212L4V

Piga (NM_011081) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Piga (NM_011081) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Piga
Synonyms:	A1194334; Pig-a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011081
ORF Size:	1458 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222212).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011081.2
RefSeq Size:	3593 bp
RefSeq ORF:	1458 bp
Locus ID:	18700
UniProt ID:	Q64323
Cytogenetics:	X 76.49 cM
Gene Summary:	Necessary for the synthesis of N-acetylglucosaminyl-phosphatidylinositol, the very early intermediate in GPI-anchor biosynthesis.[UniProtKB/Swiss-Prot Function]



[View online »](#)