

Product datasheet for **MR222222L3V**

Ctcf1 (NM_001081387) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ctcf1 (NM_001081387) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ctcf1
Synonyms:	Boris; OTTMUSG00000016680
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001081387
ORF Size:	1908 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222222).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001081387.2 , NP_001074856.1
RefSeq Size:	3217 bp
RefSeq ORF:	1911 bp
Locus ID:	664799
UniProt ID:	A2APF3
Cytogenetics:	2 H3



[View online »](#)

Gene Summary:

Testis-specific DNA binding protein responsible for insulator function, nuclear architecture and transcriptional control, which probably acts by recruiting epigenetic chromatin modifiers. Plays a key role in gene imprinting in male germline, by participating in the establishment of differential methylation at the IGF2/H19 imprinted control region (ICR). Directly binds the unmethylated H19 ICR and recruits the PRMT7 methyltransferase, leading to methylate histone H4 'Arg-3' to form H4R3sme2. This probably leads to recruit de novo DNA methyltransferases at these sites. Seems to act as tumor suppressor. In association with DNMT1 and DNMT3B, involved in activation of BAG1 gene expression by binding to its promoter. Required for dimethylation of H3 lysine 4 (H3K4me2) of MYC and BRCA1 promoters.[UniProtKB/Swiss-Prot Function]