

Product datasheet for **MR222234**

Osbp19 (NM_173350) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Osbp19 (NM_173350) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Osbp19
Synonyms:	2600011I06Rik; AU015843; Orp-9
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR222234 representing NM_173350
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGCTCTCTGTGCGCCGCTTGCAGGTTGGATTTCAGGATTCATCCCAGTGTCCAAGACTTTGATA
 AGAAACTTACCGAGGCTGACGCGTACCTGCAGATCTTGATAGAACAATTAAGCTTTTTGATGACAAGCT
 TCAAAATTGTAAAGATGATGAACAGAGAAAGAAAGTTGAAACTCTCAAAGACACAACAAATAGCATGGTA
 GAATCAATTAACACTGCATTGTGTTGCTACAGATTGCTAAAGACCAGAGTAATGCGGAGCAGCAGCGAG
 ATGGAATTATAAGTACTATTAATCCTGTAGATGCAATATACCAGCCTAGTCCCTTGGAACTGTGATCAG
 CACAATGCCTTCCAGACTGCCTTACCTCCAGAACCCGCTCAGTTGTGTAAGTCAGAGCAGCGTCCATCT
 TCCTTACCTGTTGACCTGTGTTAGCTACCTGGGACATCATCAGACTCCAACACCAAAATAGTACAGGCA
 GTGGAACTCACCACCTAGCAGCAGTCTGACTCCTCCAGCCATGTCAACTTGTCTCAAATACAGTCCC
 AGAGTTCTCTTACTCTAGCAGTGAAGATGAGTTCTATGATGCTGATGAATTCCATCAAAGTGGCTCGTCC
 CCAAAGCGCTTAATAGATTCTTCTGGATCTGCCTCAGTCTTGACACACAGCAGCTCCGGAATAGCTTAA
 AACGCCAGATACCACAGAGTCTCTGAATTCCTCCATGTCCAATGGCACAAAGCGATGCTGATCTTTTTGA
 CTCACATGACGACAGAGATGATGATGGGGAGGCTGGGTGAGTGGAGGAGCACAAGAGCGTTATCATGCAC
 CTCTTATCACAAGTCAGGCTGGGGATGGACCTCACAAAGGTAGTTCTTCCAACGTTTATTTCTCGAGAGAA
 GATCTCTGTTAGAAATGTATGCAGACTTTTTCGCACATCCAGACCTGTTCTGAGCATTAGTGTATGAGAA
 GGATCCAGGGATCGAATGGTTCAGGTTGTGAAATGGTACCTCTCGGCCCTCCATGCAGGAAGGAGAGGA
 TCGGTGGCCAAAAGCCGTACAATCCTATTTGGGTGAGATCTTTCAGTGTCACTGGACGTTGCCGAATG
 ATACTGAAGAGAACGCAGAGCTCGTTTCAGAAGGGCCGGTCCCTGGGTTTCTAAGAACAGTGAACATT
 TGTGGCTGAGCAAGTTTCCACCATCCGCCATTTTCAGCCTTTTATGCTGAGTGTTTTAAACAAGAAGATA
 CAATTCATGCTCATATCTGGACTAAATCAAAATTCCTTGGGATGTCAATTGGGGTACACAACATAGGTC
 AGGGCTGTGTCTCGTGTCTGGAGTACGATGAGCACTACATCCTCACGTTCCCAATGGCTATGGAAGGTC
 TATCCTGACAGTGCCTGGGTGGAATGGGAGGAGAAATGCAATATCAACTGCTCCAAAACGGGTTACAGC
 GCAAACATCGTCTTCCACTAAGCCTTCTATGGGGCAAGAAGCACAGAATTACTGCAGAGATTTTTT
 CTCCGAATGACAAGAAATCCTTCTGCTCAATTGAAGGGGAATGGAATGGTATCATGTATGCAAAATACGC
 AACAGGGGAAAACACTGTCTTTGTAGACACCAAGAAATGCCTATAATCAAGAAAAGGTGAGGAAGTTG
 GAAGATCAGAATGAGTATGAGTCCCGCAGCCTTTGGAAGGATGTCACTTTCAATTTAAAAATCAGAGACA
 TTGATGCAGCAACGGAAGCAAAGCACAGACTTGAAGAAAAGACAAAGAGCAGAAGCCCGAGAAAGGAAGGA
 GAAGGAAATTCAGTGGGAGACGAGGCTTTTACGAGGATGGCGAATGCTGGGTTTACGATGAGCCTTTA
 CTGAAGCGTCTTGGTGTGTGAAGCAT

AG**GCGACCG**ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR222234 representing NM_173350
 Red=Cloning site Green=Tags(s)


MALLLAACGGLDSGFIPSVQDFDKKLTEADAYLQILIEQLKLFDDKLQNCCKDDEQRKKVETLKDTTNSMV
 ESIKHCIVLLQIAKDQSNAEQHADGIIISTINPVDIYQPSLEPVISTMPSQTALPPEPAQLCKSEQRPS
 SLPVGPVLATLGHHQTPNSTGSGNSPPSSSLTPPSHVNLSPNTVPEFSYSSSEDEFYDADEFHQSGSS
 PKRLIDSSGSASVLTHSSSGNSLKRPDTTESLNSMSNGTSDADLFDSHDDRDDDGEAGSVEEHKSVMH
 LLSQVRLGMDLTKVVLPTFILERRSLEMYADFFAHPDLFVSI SDQKDPDRMVQVVKWYLSAFHAGRRG
 SVAKKPYNPILGEIFQCHWTL PNDTEENAELVSEGPVPVWSKNSVTFVAEQVSHHPPISAFYAECFNKKI
 QFNAHIWTKSKFLGMSIGVHNIQGQCVSCLYEDEHYILTFPNGYGRSILTPPWVELGGECNINCSKTGYS
 ANIVFHTKPFYGGKKHRITAEIFSPNDKKSFCISIEGEWNGIMYAKYATGENTVFVDTKKLP I IKKKVRKL
 EDQNEYESRSLWKDVTFLNKIRDIDAATEAKHRL EERQRAE ARERKEKEIQWETRLFHEDGECWVYDEPL
 LKRLGAVKH


SGP TRRRLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-RsrII

Cloning Scheme:



Plasmid Map:


ACCN: NM_173350

ORF Size: 1917 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_173350.3](#), [NP_775485.2](#)

RefSeq Size: 3148 bp
RefSeq ORF: 1920 bp
Locus ID: 100273
UniProt ID: [A2A8Z1](#)
Cytogenetics: 4 C7
MW: 72 kDa