

Product datasheet for **MR222235**

Osbp19 (NM_001134791) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Osbp19 (NM_001134791) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Osbp19
Synonyms:	2600011I06Rik; AU015843; Orp-9
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR222235 representing NM_001134791
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCTCTCTGCTGGCCGCTTGCAGAGTTTGGATTCAAGATTATCCCAGTGTCCAAGACTTTGATA
 AGAAACTTACCGAGGCTGACGCGTACCTGCAGATCTTGATAGAACAAATTAAGCTTTTTGATGACAAGCT
 TCAAAATTGTAAGATGATGAACAGAGAAAGAAAGTTGAAACTCTCAAAGACACAACAAATAGCATGGTA
 GAATCAATTAACACTGCATTGTGTTGTACAGATTGCTAAAAGTACTATTAATCCTGTAGATGCAATAT
 ACCAGCCTAGTCCCTTGGAACTGTGATCAGCACAATGCCTTCCAGACTGCCTTACCTCCAGAACCCTGC
 TCAGTTGTGTAAGTCAGAGCAGCGTCCATCTTCTTACCTGTTGGACCTGTGTTAGCTACCTTGGGACAT
 CATCAGACTCCAACCAAAATAGTACAGGCAGTGGAACTCACCACTAGCAGCAGTCTGACTCTCCCA
 GCCATGTCAACTTGTCTCCAAATACAGTCCCAGAGTCTTACTCTAGCAGTGAAGATGAGTTCTATGA
 TGCTGATGAATCCATCAAAGTGGCTCGTCCCCAAAGCGCTTAATAGATTCTTCTGGATCTGCCTCAGTC
 TTGACACACAGCAGCTCCGAAAATAGCTTAAAACGCCAGATACCACAGAGTCTCTGAATTCCTCCTATGT
 CCAATGGCACAAGCGATGCTGATCTTTTTGACTCACATGACGACAGAGATGATGATGGGGAGGCTGGGTC
 AGTGGAGGAGCACAAGAGCGTTATCATGCACCTCTTATCAAGTCAGGCTGGGGATGGACCTCACAAAAG
 GTAGTTCTTCCACGTTTATTCTCGAGAGAAGATCTCTGTTAGAAATGTATGCAGACTTTTTCGCACATC
 CAGACCTGTTTCGTGAGCATTAGTATCAGAAGGATCCCAGGGATCGAATGGTTCAGGTTGTGAAATGGTA
 CCTCTCGGCCTTCCATGCAGGAAGGAGAGGATCGGTGGCCAAAAGCCGTACAATCCTATTTTGGGTGAG
 ATCTTTTCAGTGTACTGGACGTTGCCGAATGATACTGAAGAGAACGCAGAGCTCGTTTCAGAAGGGCCGG
 TTCCCTGGGTTTCTAAGAACAGTGTAAACATTTGTGGCTGAGCAAGTTTCCACCATCCGCCATTTTCAGC
 CTTTTATGCTGAGTGTTTTAACAAGAAGATACAATTCAATGCTCATATCTGGACTAAATCAAATTCCTT
 GGGATGTCAATTGGGGTACACAACATAGGTCAGGGCTGTGTCTGTCTGGAGTACGATGAGCACTACA
 TCCTCACGTTCCCAATGGCTATGGAAGGTCTATCCTGACAGTGCCTGGGTGGAATGGGAGGAGAATG
 CAATATCAACTGCTCCAAAACGGGTTACAGCGCAAACATCGTCTTCCACACTAAGCCTTTCTATGGGGGC
 AAGAAGCACAGAATTACTGCAGAGATTTTTTCTCCGAATGACAAGAAATCCTTCTGCTCAATTGAAGGGG
 AATGGAATGGTATCATGTATGCAAAATACGCAACAGGGGAAAACACTGTCTTTGTAGACACCAAGAAAT
 GCCTATAATCAAGAAAAGGTGAGGAAGTTGGAAGATCAGAATGAGTATGAGTCCCGCAGCCTTTGGAAG
 GATGTCACTTTCAATTTAAAAATCAGAGACATTGATGCAGCAACGGAAGCAAAGCACAGACTTGAAGAAA
 GACAAAGAGCAGAAGCCCGAGAAAGGAAGGAGAAGGAAATTCAGTGGGAGACGAGGCTCTTTCACGAGGA
 TGGCGAATGCTGGGTTTACGATGAGCCTTTACTGAAGCGTCTTGGTGTGTGAAGCAT

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>MR222235 representing NM_001134791
 Red=Cloning site Green=Tags(s)

MALLLAACGGLDSGFIPSVQDFDKKLEADAYLQILIEQLKLFDDKLQNCCKDDEQRKKVETLKDTTNSMV
 ESIKHCIVLLQIAKSTINPVDATYQPSLEPVI STMPSTALPPEPAQLCKSEQRPSSLVGPVLAATLGH
 HQTPTPNSTGSGNSPPSSSLTPPSHVNLSPNTVPEFSYSSSEDEFYDADEFHQSGSSPKRLIDSSGSASV
 LTHSSGNSLKRPDTTESLNSSMSGTSDADLFDSHDDRDDGEGSVEEHKSVIMHLLSQVRLGMDLTK
 VVLPTFILERRSLEMYADFFAHPDLFVSI SDQKDRDRMVQVVKWYLSAFHAGRRGSVAKKPYNPILGE
 IFQCHWTLPNDEENAELVSEGPVPWWSKNSVTFVAEQVSHPPISAFYAEFCFNKIKQNAHIWTKSKFL
 GMSIGVHNIGQGCVCLEIDEHYILTFPNYGRSILVTPWVELGGEENINCSKTGYSANIVFHTKPFYGG
 KKHRITAEIFSPNDKKSFCSEIEGEWNGIMYAKYATGENTVFDTKKLP IIKKKVRKLEDQNEYESRSLWK
 DVTFNLKIRDIDAATEAKHRLERQRAEARERKEKEIQWETRLFHEDGECWVYDEPLLKRLGAVKH


SGPTRTRRLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_001134791

ORF Size: 1878 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001134791.2 , NP_001128263.1
RefSeq Size:	3109 bp
RefSeq ORF:	1881 bp
Locus ID:	100273
Cytogenetics:	4 C7
MW:	70.6 kDa