

Product datasheet for **MR222248**

Dynll2 (NM_001168471) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dynll2 (NM_001168471) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Dynll2
Synonyms: 1700064A15Rik; 6720463E02Rik; C87222; Dlc2; DLC8; DLC8b
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR222248 representing NM_001168471
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTGACCGGAAGGCAGTGATCAAGAACGCAGACATGTCTGAGGACATGCAACAGGATGCCGTTGACT
GTGCCACACAGGCCATGGAGAAGTACAACATAGAGAAGGACATTGCTGCCTATATCAAGAAGGAATTTGA
CAAGAAATATAACCCTACCTGGCATTGTATCGTGGCCGGAATTTGGCAGCTACGTCACACACGAGACA
AAGCACTTCATCTATTTTACTTGGTCAAGTTGCAATCCTCCTTCAAGTCGGGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR222248 representing NM_001168471
Red=Cloning site Green=Tags(s)
MSDRKAVIKNADMSEDMQQDAVDCATQAMEKYNIEKDIAAYIKKEFDKYNPTWHCIVGRNFGSYVTHET
KHFIYFYLGQVAILLFKSG


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_001168471

ORF Size: 267 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001168471.1](#), [NP_001161943.1](#)

RefSeq Size: 2270 bp

RefSeq ORF: 270 bp

Locus ID: 68097


UniProt ID: [Q9D0M5](#)

Cytogenetics: 11 C

MW: 10.3 kDa

Gene Summary: Acts as one of several non-catalytic accessory components of the cytoplasmic dynein 1 complex that are thought to be involved in linking dynein to cargos and to adapter proteins that regulate dynein function. Cytoplasmic dynein 1 acts as a motor for the intracellular retrograde motility of vesicles and organelles along microtubules. May play a role in changing or maintaining the spatial distribution of cytoskeletal structures (By similarity). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR222248