

## Product datasheet for **MR22258L4V**

### Eda2r (NM\_001161432) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Eda2r (NM_001161432) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Eda2r
Synonyms:	9430060M22Rik; TNFRSF27; Xedar
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001161432
ORF Size:	894 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR22258).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001161432.1</a> , <a href="#">NP_001154904.1</a>
RefSeq Size:	4798 bp
RefSeq ORF:	894 bp
Locus ID:	245527
UniProt ID:	<a href="#">Q8BX35</a>
Cytogenetics:	X C3
Gene Summary:	Receptor for EDA isoform A2, but not for EDA isoform A1. Mediates the activation of the NF-kappa-B and JNK pathways. Activation seems to be mediated by binding to TRAF3 and TRAF6 (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)