

Product datasheet for **MR222274L2V**

Dstn (NM_019771) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dstn (NM_019771) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dstn
Synonyms:	2610043P17Rik; ADF; AU042046; corn1; Dsn; sid23p
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_019771
ORF Size:	495 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222274).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019771.2 , NP_062745.1
RefSeq Size:	1918 bp
RefSeq ORF:	498 bp
Locus ID:	56431
UniProt ID:	Q9R0P5
Cytogenetics:	2 70.89 cM
Gene Summary:	Actin-depolymerizing protein. Severs actin filaments (F-actin) and binds to actin monomers (G-actin). Acts in a pH-independent manner.[UniProtKB/Swiss-Prot Function]



[View online »](#)