

Product datasheet for MR22281

Cbfb (NM_001161457) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cbfb (NM_001161457) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Cbfb
Synonyms: AI893578; PEA2; Pebp2; PEBP2b; Pebpb2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR22281 representing NM_001161457
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCGCGTCGTCCTCCGGACCAGAGGAGCAAGTTCGAGAACGAGGAGTCTTCAGGAAGCTGAGCCGG
AGTGCGAGATTAAGTACACGGCTTCAGGGACCGCCACGAGGAGCGCCAGACAGCTCCAGAACGC
CTGCCGACGGTCGCTCGGAGATCGCTTTTGTGGCTACAGGAACCAATCTGTCTCTCCAGTTTTTCCG
GCCAGCTGGCAGGAGAACAGCGACAAACACCTAGCCGGGAATATGTCGACTTAGAGAGAGAAGCAGGCA
AGGTATACTTGAAGGCTCCCATGATTCTGAATGGAGTGTGTTATATGGAAGGGCTGGATTGATCTCCA
CAGATTGGATGGTATGGGTTGCCTGGAGTTTATGAGGAGCGAGCCAGGCAAGAAGACAGCAAGACCCT
AGTCTGGTTCTAACTTAGGTGGCGGTGATGATCTCAAATTCGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR22281 representing NM_001161457
 Red=Cloning site Green=Tags(s)

MPRVVPDQRSKFENEFFRKL SRECEIKYTGFDRPHEERQTRFQNA CRDGRSEIAFVATGTNLSLQFFP
ASWQGEQRQTPSREYVDLEREAGKVYLKAPMILNGVCVIWKGWIDLHRLDGMGCLFDEERAQARRQDP
SPGSNLGGDDLKLR

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001161457.1, NP_001154929.1</u>
RefSeq Size:	2797 bp
RefSeq ORF:	468 bp
Locus ID:	12400
UniProt ID:	<u>Q08024</u>
Cytogenetics:	8 53.04 cM
MW:	18.5 kDa
Gene Summary:	Forms the heterodimeric complex core-binding factor (CBF) with RUNX family proteins (RUNX1, RUNX2, and RUNX3). RUNX members modulate the transcription of their target genes through recognizing the core consensus binding sequence 5'-TGTGGT-3', or very rarely, 5'-TGCGGT-3', within their regulatory regions via their runt domain, while CBF is a non-DNA-binding regulatory subunit that allosterically enhances the sequence-specific DNA-binding capacity of RUNX. The heterodimers bind to the core site of a number of enhancers and promoters, including murine leukemia virus, polyomavirus enhancer, T-cell receptor enhancers, LCK, IL3 and GM-CSF promoters (Probable). CBF complexes repress ZBTB7B transcription factor during cytotoxic (CD8+) T cell development. They bind to RUNX-binding sequence within the ZBTB7B locus acting as transcriptional silencer and allowing for cytotoxic T cell differentiation (PubMed:18258917).[UniProtKB/Swiss-Prot Function]