

Product datasheet for **MR222285**

Gab2 (NM_001162477) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gab2 (NM_001162477) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Gab2
Synonyms:	A1463667; D130058117Rik; p97
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR222285 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGCGGCGCGCGCGGACGACGTGGTGTGTACCGGCTGGCTGAGAAATCGCCTCCCGAGAAGAAGT
TGAGGCGCTATGCCTGGAAGAAACGCTGGTTTATACTTCGGAGTGGCCGAATGAGTGGAGATCCAGATGT
TCTGGAATACTACAAGAATGAGCACTCCAAGAAACCCCTGCGGATCATCAACCTGAACCTGTGTGAGCAG
GTGGATGCAGGCCTGACCTTCAACAAGAAAGAGCTGCAGGATAGTTTTGTGTTTATCAAGACCAGCG
AGCGCACATTTTACCTGGTGGCTGAGACAGAGGCTGACATGAATAAGTGGGTCCAGAGCATCTGCCAGAT
CTGCGGCTTCAATCAGGCTGAAGAGAGCACAGACTCCCTGAGGAACCTTTCTTCAGCCAGTCATGGTCCC
CGCTCTTCTCCAGCTGAGTTCAGCTCCAGTCAGCACCTGCTCCGAGAACGGAAGTCTCAGCCCCCTCAC
ACTCTAGCCAGCCTACTTTATTCAGTTTGAGCCCCCTGTGTCAAGCCACATGCAGCCTACCTTGTCCAC
CAGTGCACCTCAGGAGTATCTCTACTTGCACCAGTGCATAAGCAGAAGGACAGAAAATGCAAGGAGTGCC
AGCTTCTCTCAGGGCACCCGGCAGAAGAGTGATACAGCTGTGCAAAAACCTTGCCCAAGAGCAATGGACACT
GTATCAACGGCGTCGGAGGTCAAGTCCATGGCTTCTATAGCCTTCCCAAGCCAAGCCGACACAATACAGA
ATTCAAAGACAGTACTTATGATCTCCCACGGAGCCTGGCTTCCCATGGCCACACCAAGAGCAGCCTCACA
GGGTCTGAGACTGATAACGAGGATGTGTACACCTTCAAGATGCCAGCAACACCTGTGTGCGGAACTTG
GAGACCTCCTTGTGGACAATATGGATGTCCAACCACTCCTCTCTCAGCCTACCAGATCCCTAGAACATT
CACACTGGACAAGAACCACAATGCCATGACAGTGGCCACTCCTGGAGATTAGCCATAGCTCCCCACCC
CGGCCACCAAGCAAGTCAGGCAGAAACCTCAATGGGGCAGCATTAGCAAAAGACCTCAATCAGCG
AAAATAGCAGATCTGTAGCTGCTACTATCCCCAGGCGCAATACCCTCCCTGCAATGGACAACAGCCGACT
CCATCGAGCTTCTTCTGTGAGACCTACGAGTACCCGGCACGAGGCAAGTGGGGAAAGTGCCAGCTGGTCT
GCTGAACCTCCAGGAAAGACTGCCGTAGGTGATCAAATAGTGCCAGCTCTGATGACAACACTACGTGCCCA
TGAACCCAGGTTCTTCTACCCTGCTGGCTATGGAACGACCAGGGGACAACCTCCAGAGTGTCTACATCCC
CATGAGCCAGGACCCCATCACTTTGACCCACTTGGCTACCCGTCCACAGCCCTTCTATTACAGAGGC
CCCAGCCGAGGAAGTGAATCCAGCCACCCCGGTCAACCGAAACCTCAAGCCTGACAGAAAAGCAAAGC
CAACACCCCTTGACCTGAGAAACAACACTGTATCGATGAGCTGCCCTTCAAGTCACTGTACCAAGTC
TTGGTCCAGGATCAACAGCCACACCTTTAACTCCAGTTCTCCAGTACTGCCGTCCTATCTCCACACAG
AGCATCACCAGCACAGACTCGGGAGACAGTGGGAGAACTATGTCCCTATGCAAAACCCAGTATCTGCAT
CCCTGTTCCCAGTGGCCTAACAGCCAGCTCCAAGAAGAGTACTGGCAGTGTGGATTATCTCGCCCT
GGACTTCCAGCCGGGCTCCCCAAGCCCTCACCGCAAGCCATCCACATCATCTGTACATCAGATGAGAAG
GTAGACTATGTCCAAGTGGATAAAGAGAAGACCCAGGCCCTGCAGAACCCATGCAGGAGTGGACAGATG
TGCGGCAGTCTCCGAACCTTCAAGGGTGCCAAGCTG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR222285 protein sequence
Red=Cloning site Green=Tags(s)


MSGGGDDVVCTGWLKSPPEKKLRRYAWKKRWFILRSGRMSGDPDVLEYKNEHSKKPLRIINLNLCEQ
VDAGLTFNKKELQDSFVFDIKTSERTFYLVAETEADMNKWVQSIQICGFNQAEESTDSLRLSSASHGP
RSSPAEFSSSQHLLRERKSSAPSHSSQPTLFTFEPPVSSHMQPTLSTAPQEYL YLHQICISRRTENARSA
SFSQGTRQKSDTAVQKLAQSNHCINGVGGQVHGFYSLPKPSRHNTFEKDYDLPRSLASHGHTKSSLT
GSETDNEVDVYTFKMPNTLCRELDLLVDNMDVPTTPLSAYQIPRTFTLDKNHNAMTVATPGDSAIAPPP
RPPKPSQAETPQWGSIQQRPPISENSRVAATIPRRNTLPAMDNSRLHRASSCETYEYPARGSGESASWS
AEPGKTAVGRSNSASSDDNYVPMNPGSSTLLAMERPGDNSQSVYIPMSPGPHFDPLGYPSTALPIHRG
PSRGSEIQPPPVRNLKPDRAKAPTPDLRNNTVIDELPFKSPVTKSWRINSHTFNSSSQYCRPISTQ
SITSTDSGDSEENYVPMQNPVSASPVSGTNSPAPKKSTGSDYLALDFQPGSPSPHRKPSTSSVTSDEK
VDYVQVDEKKTQALQNTMQEWDVQRSSSEPSKGAKL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001162477

ORF Size: 2001 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001162477.1](#), [NP_001155949.1](#)

RefSeq Size: 6005 bp


RefSeq ORF: 1998 bp

Locus ID: 14389

Cytogenetics: 7 E1

MW: 73.3 kDa

Product images:



Circular map for MR222285