

## Product datasheet for **MR222289L4V**

### Dap3 (NM\_001164533) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Dap3 (NM_001164533) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dap3
Synonyms:	4921514D13Rik; DAP-3; MRP-S29; S29mt
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001164533
ORF Size:	1098 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222289).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001164533.1</a> , <a href="#">NP_001158005.1</a>
RefSeq Size:	4538 bp
RefSeq ORF:	1101 bp
Locus ID:	65111
UniProt ID:	<a href="#">Q9ER88</a>
Cytogenetics:	3 39.01 cM
Gene Summary:	Involved in mediating interferon-gamma-induced cell death.[UniProtKB/Swiss-Prot Function]



[View online »](#)