

Product datasheet for **MR222302L4V**

Cysltr2 (NM_001162412) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cysltr2 (NM_001162412) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cysltr2
Synonyms:	2300001H05Rik; Cltr2; Cysl2; CYSLT2R
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001162412
ORF Size:	927 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222302).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001162412.1 , NP_001155884.1
RefSeq Size:	1745 bp
RefSeq ORF:	930 bp
Locus ID:	70086
UniProt ID:	Q920A1
Cytogenetics:	14 D3
Gene Summary:	Receptor for cysteinyl leukotrienes. The response is mediated via a G-protein that activates a phosphatidylinositol-calcium second messenger system. The rank order of affinities for the leukotrienes is LTC4 = LTD4 >> LTE4.[UniProtKB/Swiss-Prot Function]



[View online »](#)